TOTAL FACE MASK FOR NIV THERAPY

HAROL Your breathing Our mission



INFORMATION SHEET Total Face Mask for NIV therapy



The **FACE-HAR** mask is used in non-invasive ventilation (CPAP / BI-LEVEL-PS) in combination with **double-limb ventilation circuits**. It is equipped with a 4-point fabric mask **headgear** and has a specifically designed **anatomical shape** to allow maximum comfort and perfect adherence to any kind of face shape.

It is equipped with a **sealed access** for nasogastric tube and a **connection fitting to the ventilation circuit** in blue color (not ventilated). An optional **fitting with anti-suffocation valve** (EE) is also available.

IT IS RECOMMENDED FOR:

Non-cardiogenic pulmonary edema (A.L.I. - ARDS) - Chest injury - Atelectasis - Post-operative hypoxemia - Asthma - Exacerbated COPD - Pre and post-extubation respiratory assistance.

>> TECHNICAL FEATURES

- The mask is **fully transparent**, and it is equipped with a **very soft** silicone lip cushion, without metal parts for MR use.
- It is available in **3 adult sizes**, to perfectly fit the patient's face.
- It is equipped with a **leak-free sealed access port** for **nasogastric and fiberscope tubes**, to avoid leaks and interruptions in ventilation.

THIS IS A DISPOSABLE AND NON-STERILE DEVICE


BENEFITS

CONVENIENCE: the **dead space is reduced** for better synchronization with the ventilator (from 300 to 530 ml, depending on the size)

PERFECT SEAL: The silicone cushion fits perfectly and firmly to the entire face, ensuring a good seal with less pressure on the patient.

SAFE MATERIAL: The device is **LATEX-FREE** and **DEHP-FREE** (it does not contain phthalates).

QUALITY: Excellent performance at an affordable price.



SPECIFICATIONS		
CODE	SIZE	PACKAGING
CPAP9013811	Adults S	Single packed
CPAP9013812	Adults M	Single packed
CPAP9013813	Adults L	Single packed

