



# MiniMe<sup>®</sup> 2

featuring AIR<sup>°</sup>gel<sup>®</sup> with Advanced Cushion Technology

Designed for children ages 2-12, our MiniMe 2™ mask features our patented AIR<sup>°</sup>gel with Advanced Cushion Technology. With increased moldability, MiniMe 2 provides your smallest patients the ultimate fit and maximum comfort.

## The Next Generation of Pediatric Masks

### AIR<sup>°</sup>gel

Maximum comfort, minimum weight. Our patented AIR<sup>°</sup>gel makes the MiniMe 2 the most gentle fitting mask on the market.



### Touchless Spacebar™

Unique design eliminates any contact with the patient's forehead.

### Size Indicator

Both Vented and Non-Vented come in two sizes. **S/M** **M/L**

### Custom Fit Technology

The shapeable shell mask can be easily molded to fit individual needs.

*Flex, Fit, Forget.*

### Advanced Cushion Technology (A.C.T.)

Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.

### 360° Rotating Elbow

For optimal tube positioning and easier mobility

### Vented and Non-Vented Versions Available

For CPAP and Bi-level therapy in children



Made in the USA

## What makes our MiniMe 2 Mask different?

- Designed for ages 2 – 12, the MiniMe 2 incorporates the latest technology including our Touchless Spacebar and Advanced Cushion Technology (A.C.T.).
- *Custom Fit Technology* - a unique soft-shell mask with a shapeable wire allowing for optimal flexibility, fit, and comfort.
- Available in two sizes with 4-point simple and breathable headgear.

## Advanced Cushion Technology (A.C.T.)

Our patented AIR<sup>o</sup>gel with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal.



## Custom Fit Technology

BEND IT.



MOLD IT.



SHAPE IT.



## MiniMe 2 Masks with Headgear Specifications:

Product Number	Style	Headgear size	Leak Port	Mask Size	Mask Weight	Dead Space	Mask Dimensions
60217	Vented	Small	✓	S/M	31.6 g	22 ml	2.32"D x 2.00"W x 1.71"H
60218	Vented	Large	✓	M/L	40.3 g	40 ml	2.59"D x 2.69"W x 2.14"H
60257	Non-Vented	Small	✗	S/M	31.6 g	22 ml	2.32"D x 2.00"W x 1.71"H
60258	Non-Vented	Large	✗	M/L	40.3 g	40 ml	2.59"D x 2.69"W x 2.14"H
55048	Headgear	Small	N/A	N/A	N/A	N/A	N/A
55051	Headgear	Large	N/A	N/A	N/A	N/A	N/A



Does not contain natural rubber latex

Also from Sleepnet:

## Veraseal 2

A Giant Step Forward in Patient Comfort



Contact us today.

To place your order, call an authorized distributor.

Call 1.800.742.3646

Visit [www.sleepnetmasks.com](http://www.sleepnetmasks.com) to learn more.

Sleepnet Corporation, 5 Merrill Industrial Drive, Hampton, New Hampshire 03842 USA



Scan code to view MiniMe 2 Fitting Tips video.