

# ECG ELECTRODE



## TELECTRODE

Connect with Confidence

Produces ECG electrodes called 'Telectrodes'  
which are used for various ECG tests

These are excellent components for sensitive monitoring

**Bio Protech Inc.**

# The perfect companion for heart health.



With our innovative solutions and products, we can assist healthcare facilities in optimizing their operational costs without compromising on the quality of patient care diagnostics.

## Best price, Best performance

Top-notch performance at a low price. An economical choice, but no compromise on performance. With help from Bioprotech, you can meet the challenge by standardizing electrodes without compromising on product performance

## Comfortable fit

Emphasizing the design and material selection of the product, the ECG electrodes can touch the skin smoothly and be worn comfortably.

## Accurate and reliable measurement

High-performance sensor and ECG electrodes of proven quality provide accurate ECG measurement. Quick and consistently excellent signal quality is obtained.

## Backing choice


You can choose from 3 types ( Foam, Cloth, Clear tape) of backing for conformability, patient comfort and low profile.



### PE Foam

The foam material provides a soft and cushioning effect that gently touches the skin. It prevents unnecessary moisture build-up on the skin and supports skin health.

- T715
- T716
- T717
- T718
- T712T
- T711T
- T711CT
- T724T
- T724CT



### Cloth

The cloth material adheres reliably to the skin and remains in place during movement. The electrodes provide the long term comfort and durability required by active patients.

- T815
- T815C
- T816
- T832
- T832C
- T834
- T834C



### Clear Tape

The clear tape material is transparent, and when the ECG electrodes are attached to the skin, the condition under the skin can be visually checked. This can be helpful in monitoring skin conditions or monitoring skin changes.

- T916
- T916C

By choosing Bio Protech as your partner, you can benefit from our expertise, cost-effective solutions, and commitment to high-quality patient care diagnostics. Let us help you navigate the challenges of reducing operating costs while maintaining excellence in healthcare.

# TELECTRODE


## Monitoring ECG Electrodes

### Multiple sizes and shapes

We offer ECG electrodes for adult patients in a variety of sizes and shapes to suit different body shapes and sizes.

### ECG Electrodes for Universal

Model	Backing	Size(mm)	Shape
T715	PE Foam	Dia. 55	Round
T716A	PE Foam	43 x 45	Teardrop
T717	PE Foam	Dia. 43	Round
T711T	PE Foam	37 x 42.5	Teardrop
T712T	PE Foam	45 x 49	Teardrop
T724T	PE Foam	36 x 42	Square
T816	Cloth	43 x 45	Teardrop
T916	Clear Tape	43 x 45	Teardrop



T716A



T724T



T916

### Small size option

The T718 electrode is compact size and design to suit the pediatric's anatomy.

### ECG Electrodes for Pediatric

Model	Backing	Size(mm)	Shape
T718	PE Foam	32 x 38	Teardrop



T718

Suggested Applications	T715	T716A	T717	T718	T711T	T712T	T724T	T816	T916
EMS/ED	○	○	○	○				○	○
ICU/CCU	○	○	○	○				○	○
Telemetry/Holter/Event					○	○	○	○	○
OR	○	○	○	○	○	○	○		
Cath Lab									
Stress Test							○		
Cardiac Rehab	○	○	○	○	○	○	○	○	○
Echo-Cardio									
Universal Foam Monitoring	○	○	○	○			○	○	○
Pediatric				○					
Solutions of Fragile skin									
Diaphoretic solutions									
Fluid resistant solutions	○	○	○	○		○	○		